





チタン製品カタログ



鉗子		
1	外科用ピンセット	T1
2	ピンセット（ドベキー血管,ジェラルド組織,ロートン硬膜,アドソン,クッシング）	T2～T5
3	ピンセット（マイクロ,グレーフェ,血管,レザノ血管,アドソン・ブラウン組織,アドソン組織 ニュージェント）	T6～T7
4	ピンセット（アダム,外科用,マイクロ結紮,デニスマイクロ,マイクロリング）	T8～T10
5	ジュエラーピンセット	T11～T12
6	マイクロピンセット（ヤサーギル）	T13～T14
7	クランプ鉗子（カーペンティア乳頭筋、オショーネシー）	T14～T15

止血鉗子（モスキート）	T17～T18
-------------	---------

血管クランプ		
1	ディートリックブルドッグ鉗子	T19
2	クランプ鉗子（デニス多目的,ベック乳児大動脈血管,スプラマニアン幼児大動脈血管 ドベキー血管クーリー血管,グローバー）	T19～T22
3	サテンスキー血管クランプ鉗子	T23
4	新生児血管鉗子（カスターネダ,ドベキー,キャッツミラー）	T24～T25

持針器		
1	マイクロ持針器（ヤコブソン,ヤコブソンヘビータイプ,カストロピエホ,ヤコブソン）	T27～T33
2	マイクロ持針器（ヤサーギル,ロートン）	T34
3	ヘガール型持針器（ライダー,フェル,クライル・ウッド,胸骨針用）	T35～T38
4	ドベキー・ダブルアクション持針器	T39
5	オルセン・ヘガール剪刀付持針器	T40

剪刀		
1	外科剪刀	T41
2	リスター包帯剪刀	T41
3	剪刀（ミニヴァンナス,ヴァンナス,ノイス,ウエストコット,ヤコブソンマイクロ,マイクロ）	T42～T44
4	マイクロ剪刀（ポッツ,ヤコブソン,ヤコブソン・ヴァンナス）	T44～T47
5	マイクロ剪刀（ヤサーギル,マリス・バイオネット,バイオネット）	T48～T50
6	レザノ剪刀 三連関節式	T51
7	外科剪刀	T52
8	メーヨー剪刀	T53

剥離子		
1	ヤサーギルマイクロ結紮ガイド,ヤサーギル・マイクロ血管メス,ヤサーギル・マイクロ剥離子 ヤサーギルマイクロディセクター,ヤコブソン・マイクロプローブ,ヤコブソン マイクロ結紮フォーク,クライエンビュール血管・神経鉤	T55～T56
2	マイクロフック,キュレット	T56
3	ヤサーギルマイクロくも膜メス,マイクロ鋭匙,ヤサーギルマイクロ剥離子 バイオネット（マイクロ血管メス,ディセクター,エレベーター,神経鉤）	T56～T58

キュレット,血管鉗子		
1	バイオネットマイクロリングキュレット	T59
2	ギャレット血管拡張器	T59
3	カルボンティ血管鉤	T60

レトラクター		
1	I.M.A.用心外膜レトラクター,ウッド両頭鉤,クーリー心房鉤（僧房弁用鉤） 新生児用胸骨用開創器,フィノシット胸骨用開創器 新生児	T61～T63
2	胸骨開胸器,フィノシット肋骨開胸器,モース胸骨開胸器,デュポスト開胸器 カルボンティ開胸器, I.M.A開胸器 内胸動脈用	T64～T68

外科用ピンセット

Without teeth

1813102 13cm

1813103 15cm

1813104 16cm

1813105 18cm

1813106 20cm

1813107 23cm

1813108 25cm



1813102



1813104



1813106


ドバーキー血管ピンセット

Flat handle, 2.0mm atraumatic tips

1801270 Straight, 15cm

1801271 Straight, 20cm

1801272 Straight, 24cm



1801271


ドバーキー血管ピンセット

Flat handle, 1.5mm atraumatic tips

1801274 Straight, 15cm

1801275 Straight, 20cm

1801276 Straight, 24cm



1801274



1801275


ドバーキー血管ピンセット

Flat handle, 1.0mm atraumatic tips

1801280 Straight, 15cm

1801281 Straight, 20cm

1801282 Straight, 24cm



1801281


ドバーキー血管ピンセット

Round handle, 2.0mm atraumatic tips

1801333 Straight, 15cm

1801335 Straight, 20cm


1801337 Straight, 24cm




1801335

ドバーキー血管ピンセット
Round handle

- 1801342** Straight, 1.0mm atraumatic tips, 15cm
- 1801343** Straight, 1.0mm atraumatic tips, 20cm
- 1801344** Straight, 1.0mm atraumatic tips, 24cm
- 1801345** Angled, 1.2mm atraumatic tips, 21cm




1801343



1801345

ドバーキー血管ピンセット
Round handle, Straight, 1.0mm atraumatic tips

- 1801350** Straight, 18cm
- 1801351** Straight, 21cm
- 1801352** Straight, 23cm




1801351



ドバーキー血管ピンセット
2.0mm atraumatic tips, Tungsten carbide coated behind
the 6mm length atraumatic jaw for secure needle pulling action


- 1801287** Straight, Flat handle, 20cm
- 1801346** Straight, Round handle, 20cm

1801287



1801346

ドバーキー血管ピンセット
1801283 Straight, 1.5mm atraumatic tips, 20cm
1801284 Straight, 2.0mm atraumatic tips, 20cm




1801284

ジェラルド組織ピンセット

1801285 Straight, 1.0mm atraumatic tips, 18cm

1801286 Straight, 1.0mm atraumatic tips, 20cm



1801286

ジェラルド組織ピンセット

1801002 Straight, 1.0mm Platforms, 18cm

1801004 Straight, 1.0mm Platforms, 20cm




1801004

ジェラルド組織ピンセット

Tungsten carbide coated tips, 1.0mm ring tips

1801241 Straight, 18cm

1801242 Straight, 20cm



1801242

組織ピンセット

1.2mm beak nose tip, with 10mm tungsten carbide coated platforms

1801347 20cm

1801348 23cm




1801348

コートン硬膜ピンセット

Tungsten carbide coated tips, 17.8cm


1801349 1.0mm tip width, 1×2teeth



組織ピンセット


1801354 0.4mm, 1×2teeth, 6mm platforms, 18cm

1801355 0.4mm, 1×2teeth, 7mm platforms, 21cm



アドソンピンセット

1801435 0.4mm, 1×2teeth, 15cm



クッシングピンセット

1801288 1mm serrated tips, Straight, Scraper end, 18cm

1801289 1.4mm serrated tips, Angled, 18cm

1801290 1mm serrated tips, Straight, 21cm



組織ピンセット

1801291 2x3teeth, 20cm



マイクロピンセット

- 1801293** 0.8mm tips, 1 x 2teeth,
Tungsten carbide coated tips, 11.5cm
- 1801294** 0.8mm serrated tips, 11.5cm

1801293



1801294



グレーフェピンセット

Tungsten carbide coated tips, Straight, 10cm

- 1801295** 0.5mm delicate tips, 1×2 teeth
- 1801296** 0.8mm tips, 1×2 teeth
- 1801297** 0.5mm delicate serrated tips
- 1801298** 0.8mm serrated tips

1801295



グレーフェピンセット

Tungsten carbide coated tips, Curved, 10cm

- 1801299** 0.5mm delicate tips, 1×2 teeth
- 1801300** 0.8mm tips, 1×2 teeth
- 1801301** 0.5mm delicate serrated tips
- 1801302** 0.8mm serrated tips

1801302



グレーフェピンセット

Tungsten carbide coated tips, Strong curved, 10cm

- 1801303** 0.8mm tips, 1×2 teeth
- 1801304** 0.5mm delicate serrated tips
- 1801305** 0.8mm serrated tips

1801303



グレーフェピンセット

Tungsten carbide coated tips, Angled, 10cm

- 1801306** 0.8mm tips, 1×2 teeth
- 1801307** 0.5mm delicate serrated tips
- 1801308** 0.8mm serrated tips

1801306



血管ピンセット

Narrow double row of fine serrations on each jaw

1801309 18cm

1801310 23cm

1801311 32cm

1801309



レザノ血管ピンセット

1801312 With 3 articulations, 28cm



アドソン・ブラウン組織ピンセット

1801313 1.8mm tips, 7×7 teeth, 12.1cm



アドソン組織ピンセット

Tungsten carbide coated tips, 12.1cm

1801315 1.0mm tip width, 1×2 teeth

1801316 1.0mm tip width, 1×2 teeth,
with tying platforms

1801317 1.3mm tip width,
2×3 teeth

1801315



アドソン組織ピンセット

Tungsten carbide coated tips, 12.1cm

1801003 1.0mm tips, without tying platforms

1801088 1.0mm tips, with tying platforms

1801003



ニュージェントピンセット

1801005 1.2mm tips, Tungsten carbide coated
tips, Angled, 10.8cm




ピンセット

Tungsten carbide coated tips

1801006	1mm tips, 16cm	1801009	2mm tips, 16cm
1801007	1mm tips, 20cm	1801010	2mm tips, 20cm
1801008	1mm tips, 24cm	1801011	2mm tips, 24cm


1801010



アダムピンセット

Very delicate jaws, with 10mm tungsten carbide coated platforms, Round handle, Straight

1801089 1.2mm tips, 20cm



外科用ピンセット

1801090 2mm tungsten carbide coated tips, 22.5cm




マイクロ結紮ピンセット

Tungsten carbide coated platforms, Round handle, 15cm

1801091 Straight, 0.2mm tips

1801092 Straight, 0.1mm tips



1801091




マイクロ結紮ピンセット

Tungsten carbide coated platforms, Round handle, 15cm

1801093 Straight, 0.3mm tips

1801094 Straight, 0.2mm tips



1801093





マイクロ結紮ピンセット

Tungsten carbide coated platforms, Round handle,15cm

1801095 Curved, 0.3mm tips

1801096 Curved, 0.2mm tips



1801095

マイクロ結紮ピンセット


Tungsten carbide coated platforms, Straight, Round handle,18cm

1801097 0.7mm tips

1801098 0.5mm tips

1801099 0.3mm tips

1801100 0.2mm tips



1801099

マイクロ結紮ピンセット


Tungsten carbide coated platforms , Curved, Round handle,18cm

1801101 0.7mm tips

1801102 0.5mm tips

1801103 0.3mm tips

1801104 0.2mm tips



1801103

マイクロ結紮ピンセット

Tungsten carbide coated platforms, Round handle

1801105 Straight 0.5mm tips, 20cm

1801106 Straight 0.3mm tips, 20cm

1801109 Straight 0.7mm tips, 23cm

1801110 Straight 0.3mm tips, 23cm

1801113 Straight 1.0mm tips, 25cm

1801114 Straight 0.5mm tips, 25cm

1801117 Straight 1.0mm tips, 30cm

1801119 Straight 1.0mm tips, 33cm

1801107 Curved 0.5mm tips, 20cm

1801108 Curved 0.3mm tips, 20cm

1801111 Curved 0.7mm tips, 23cm

1801112 Curved 0.3mm tips, 23cm

1801115 Curved 1.0mm tips, 25cm

1801116 Curved 0.5mm tips, 25cm

1801118 Curved 1.0mm tips, 30cm

1801120 Curved 1.0mm tips, 33cm



1801105



1801107

マイクロ結紮ピンセット

With 10mm tungsten carbide coated platforms, Straight

1801121 0.2mm tips, 20cm

1801122 0.3mm tips, 23.5cm

1801123 0.3mm tips, 27.5cm

1801123



デニスマイクロピンセット

Tungsten carbide coated platforms

1801133 Straight, 0.3mm tips, 15cm

1801138 Straight, 0.5mm tips, 18cm

1801139 Straight, 0.3mm tips, 18cm

1801145 Straight, 0.5mm tips, 21cm

1801146 Straight, 0.3mm tips, 21cm

1801135 Curved, 0.3mm tips, 15cm

1801142 Curved, 0.5mm tips, 18cm

1801143 Curved, 0.3mm tips, 18cm

1801147 Curved, 0.5mm tips, 21cm

1801148 Curved, 0.3mm tips, 21cm

1801139



1801143



マイクロリングピンセット

Tungsten carbide coated tips, 1.0mm ring tips, Round handle

1801210 15cm

1801214 18.1cm

1801215 20cm


1801216 21cm

1801217 23cm

1801218 25cm

1801219 27.5cm

1801210



1801219



マイクロリングピンセット


Tungsten carbide coated tips, Round handle

1801220 1.3mm ring tips, 15cm

1801221 1.3mm ring tips, 18cm

1801228 1.5mm ring tips, 20cm

1801230 1.5mm ring tips, 23cm



1801228



ジュエラーピンセット

1801012 1# Straight, 0.2×0.12mm tips,
12cm



1801013 1# Straight, 0.1×0.06mm tips,
12cm



1801014 2# Straight, 0.34×0.14mm tips,
12cm



1801015 2# Straight, 0.17×0.07mm tips,
12cm



1801016 3# Straight, 0.17×0.1mm tips,
12cm



1801017 3# Straight, 0.08×0.04mm tips,
12cm



1801018 3c# Straight, 0.08×0.04mm
tips, 11cm



1801019 4# Straight, 0.13×0.08mm
tips, 11cm



ジュエラーピンセット

1801020 4# Straight, 0.06×0.02mm tips,
11cm



1801022 5# Straight, 0.05×0.01mm tips,
11cm



1801023 5/45°# Bend, 0.05×0.01mm
tips, 11cm



1801025 7# Curved, 0.17×0.1mm tips,
11.7cm



1801026 0.15×0.12mm tips, 15.5cm



1801027 5a#, 0.1×0.06mm tips, 11.7cm



1801028 5a#, 0.13×0.08mm tips, 11.7cm




ヤサーギルマイクロピンセット

Tungsten carbide coated tips, Bayonet style, Straight

1801600 0.6mm tips, 21.6cm

1801604 0.9mm tips, 24cm

1801606 0.9mm tips, 25cm



1801600

ヤサーギルマイクロピンセット

1801602 Straight, 1.0mm tips, 23cm


1801602A Straight, 0.5mm tips, 23cm

1801650 Curved, 1.0mm tips, 23cm

1801650A Curved, 0.5mm tips, 23cm

1801651 Angle up, 1.0mm tips, 23cm

1801651A Angle up, 0.5mm tips, 23cm



1801650




1801651

マイクロピンセット

Bayonet style, Straight, Tungsten carbide coated tips



1801608 0.4mm tips, 26.3cm



ヤサーギルマイクロピンセット

Tungsten carbide coated tips

- 1801610** 1.3mm tips,
1×2teeth,
21.6cm
- 1801610A** 0.7mm tips,
1×2teeth, 21.6cm



1801610

ヤサーギルマイクロピンセット


Bayonet style, Straight, Tungsten carbide coated tips

- 1801620** 1.0mm ring tips, 23cm



カーペンティア乳頭筋クランプ鉗子


- 1800452** 17cm
- 1800440** 24cm
- 1800441** 27cm



1800441

オシヨーネシークランプ鉗子

Ring handle



Right Angle

1800435B 18.5cm
1800436B 20cm
1800437B 22cm
1800438B 24cm
1800439B 28cm

1/2 Curved

1800435A 18.5cm
1800436A 20cm
1800437A 22cm
1800438A 24cm
1800439A 28cm

Curved

1800435 18.5cm
1800436 20cm
1800437 22cm
1800438 24cm
1800439 28cm

モスキート止血鉗子

Very delicate serrated jaws

1814115 0.6mm tips,
Straight,12.5cm




1814116 0.6mm tips,
Slightly curved,12.5cm



1814117 0.6mm tips,
Curved,12.5cm



1814118 0.6mm tips,
Right angle,12.5cm




1814315 Straight, 0.6mm tip width,
1×2teeth,12.5cm




1814316 Slightly curved, 0.6mm tip width,
1×2teeth,12.5cm



1814317 Curved, 0.6mm tip width,
1×2teeth,12.5cm



1814318 Right angle, 0.6mm tip width
1×2teeth,12.5cm



モスキート止血鉗子

1800404 Straight, Serrated jaws, 12.5cm




1800405 Curved, Serrated jaws, 12.5cm



1800406 Straight, 1×2teeth, 12.5cm



1800407 Curved, 1×2teeth, 12.5cm



1800405

止血鉗子 Serrated jaws

1814113 Straight, 14cm




1814133 Curved, 14cm



1814313 Straight, 1×2teeth, 14cm



1814333 Curved, 1×2teeth, 14cm



1814333

ディートリックブルドック鉗子

Cross-action, Atraumatic tips, Tension 180gms

1820450 Straight, 2×11mm jaws, 55mm

1820455 Straight, 2×15mm jaws, 61mm

1820470 Straight, 2×20mm jaws, 69mm

1821450 Angled, 2×11mm jaws, 51mm

1821455 Angled, 2×15mm jaws, 57mm

1821470 Angled, 2×20mm jaws, 65mm

1821450



1821455



ディートリックブルドック鉗子

Cross-action, Serrated tips, Tension 80gms

1822450 Angled, 2×11mm jaws, 51mm

1822455 Angled, 2×15mm jaws, 57mm

1822470 Angled, 2×20mm jaws, 65mm

1823450 Straight, 2×11mm jaws, 55mm

1823455 Straight, 2×15mm jaws, 61mm

1823470 Straight, 2×20mm jaws, 69mm

1822450



1822455



1823455




デニス多目的クランプ鉗子

1800280 6cm atraumatic jaw, 18.5cm

1800288 8cm atraumatic jaw, 20cm

1800280




1800280

1800288


ベック乳児大動脈血管クランプ鉗子

1800366 Miniature aorta clamp, Atraumatic jaw, 29mm jaw length,
7mm jaw depth, 16.5cm




スプラマニアン幼児大動脈血管クランプ鉗子

1800368 Miniature aorta clamp, Atraumatic jaw, 30mm jaw length,
7mm jaw depth, 16cm



ドバーキ血管鉗子


1800369 76mm jaw length, 14mm jaw depth, 18.8cm



ドバーキ血管鉗子


1800370 Kay aorta exclusion clamp, Atraumatic jaw, 51mm jaw length, 32mm jaw depth, 20cm

1800372 Lambert-Kay aorta exclusion clamp, Atraumaic jaw, 51mm jaw length, 20mm jaw depth, 20cm




クーリー血管鉗子


1800374 Atraumatic jaw, 27cm




1800376 ドベーキ血管鉗子 Atraumatic jaw, 27cm



1800378 グローバー鉗子 23.2cm




1800388 ドベーキ血管鉗子 27cm



サテンスキー血管鉗子

- 1800380 Satinsky Clamp**
30mm jaw, 23.2cm
- 1800382 Satinsky Clamp**
24mm jaw, 23.2cm
- 1800384 Satinsky Clamp**
18mm jaw, 23.2cm



1800380

1800382








1800384

1800382




カスタネダ新生児血管鉗子

- 1800442** Straight jaw, 13.5cm
- 1800443** 30° angled jaw, 13.5cm
- 1800444** 50° angled jaw, 13.5cm
- 1800445** 70° angled jaw, 12cm
- 1800446** 90° angled jaw, 11cm
- 1800447** Curved jaw, 12cm
- 1800448** Curved jaw, 13.5cm
- 1800449** Curved jaw, 15cm
- 1800450** Curved jaw, 15cm
- 1800451** Curved jaw, 15cm



ドバーキ新生児血管鉗子

- 1800610** Straight shanks, Straight jaw, 11.5cm
- 1800611** Slightly angled shanks, Straight jaw, 11.5cm
- 1800614** Straight shanks, 60° angled jaw, 10.5cm
- 1800616** Slightly angled shanks, 6mm jaw, 10.5cm
- 1800618** Straight shanks, 6mm spoon jaw, 10.5cm
- 1800620** Slightly angled shanks, 9mm jaw, 10.5cm



1800620



1800618

キャッツミラー新生児血管鉗子

- 1800628** Straight shanks, Left angled jaw, 23mm jaw length, 11cm
- 1800630** Straight shanks, Right angled jaw, 23mm jaw length, 11cm



1800628



1800628

1800630

マイクロ持針器

Round handle, Tungsten carbide coated tips, Straight, 13.8cm

- 1802031** Without lock, 0.3mm tips
- 1802035** Without lock, 0.5mm tips
- 1802039** Without lock, 0.7mm tips
- 1802043** Without lock, 1.0mm tips
- 1802047** Without lock, 1.5mm tips

- 1802032** With lock, 0.3mm tips
- 1802036** With lock, 0.5mm tips
- 1802040** With lock, 0.7mm tips
- 1802044** With lock, 1.0mm tips
- 1802048** With lock, 1.5mm tips



1802031

マイクロ持針器

Round handle, Tungsten carbide coated tips, Curved, 13.8cm

- 1802033** Without lock, 0.3mm tips
- 1802037** Without lock, 0.5mm tips
- 1802041** Without lock, 0.7mm tips
- 1802045** Without lock, 1.0mm tips
- 1802049** Without lock, 1.5mm tips

- 1802034** With lock, 0.3mm tips
- 1802038** With lock, 0.5mm tips
- 1802042** With lock, 0.7mm tips
- 1802046** With lock, 1.0mm tips
- 1802050** With lock, 1.5mm tips



1802033

マイクロ持針器

Round handle, Tungsten carbide coated tips, Straight, 18.5cm,

- 1802092M** 0.3mm tips, Multilock
- 1802094M** 0.5mm tips, Multilock
- 1802096M** 0.7mm tips, Multilock
- 1802098M** 1.0mm tips, Multilock
- 1802100M** 1.5mm tips, Multilock

- 1802092L** 0.3mm tips, with lock
- 1802094L** 0.5mm tips, with lock
- 1802096L** 0.7mm tips, with lock
- 1802098L** 1.0mm tips, with lock
- 1802100L** 1.5mm tips, with lock



1802094M



1802098L

マイクロ持針器

Round handle, Tungsten carbide coated tips, Straight, Multilock

1802132M	0.3mm tips, 20cm
1802134M	0.5mm tips, 20cm
1802136M	0.7mm tips, 20cm
1802138M	1.0mm tips, 20cm
1802140M	1.5mm tips, 20cm
1802194M	0.5mm tips, 23cm
1802196M	0.7mm tips, 23cm
1802198M	1.1mm tips, 23cm
1802200M	1.5mm tips, 23cm
1802214M	0.5mm tips, 25cm
1802216M	0.7mm tips, 25cm
1802218M	1.1mm tips, 25cm
1802220M	1.5mm tips, 25cm
1802234M	0.5mm tips, 30cm
1802236M	0.7mm tips, 30cm
1802238M	1.1mm tips, 30cm
1802240M	1.5mm tips, 30cm

1802136M



1802198M



マイクロ持針器

Round handle, Tungsten carbide coated tips, Curved, Multilock

1802102M 0.3mm tips, 18.5cm
1802104M 0.5mm tips, 18.5cm
1802106M 0.7mm tips, 18.5cm
1802108M 1.0mm tips, 18.5cm
1802110M 1.5mm tips, 18.5cm

1802142M 0.3mm tips, 20cm
1802144M 0.5mm tips, 20cm
1802146M 0.7mm tips, 20cm
1802148M 1.0mm tips, 20cm
1802150M 1.5mm tips, 20cm

1802204M 0.5mm tips, 23cm
1802206M 0.7mm tips, 23cm
1802208M 1.1mm tips, 23cm
1802210M 1.5mm tips, 23cm

1802224M 0.5mm tips, 25cm
1802226M 0.7mm tips, 25cm
1802228M 1.1mm tips, 25cm
1802230M 1.5mm tips, 25cm

1802244M 0.5mm tips, 30cm
1802246M 0.7mm tips, 30cm
1802248M 1.1mm tips, 30cm
1802250M 1.5mm tips, 30cm

1802104M



1802146M



マイクロ持針器

Round handle, Tungsten carbide coated tips, Multilock

1802286M 0.7mm tips, Straight, 27.5cm

1802296M 0.7mm tips, Curved, 27.5cm



ヤコブソン持針器

Round handle, Tungsten carbide coated tips, Straight, for large needles, Multilock

1802300M 1.5 × 17.5mm jaw, 17.8cm

1802301M 1.5 × 17.5mm jaw, 19cm

1802302M 1.5 × 17.5mm jaw, 20.5cm

1802303M 1.5 × 17.5mm jaw, 22cm

1802304M 1.5 × 17.5mm jaw, 23cm

1802305M 1.5 × 17.5mm jaw, 25cm



ヤコブソン持針器

Round handle, Tungsten carbide coated tips, Curved, for large needles, Multilock

1802306M 1.5 × 17.5mm jaw, 17.8cm

1802308M 1.5 × 17.5mm jaw, 20.5cm

1802309M 1.5 × 17.5mm jaw, 22cm

1802311M 1.5 × 17.5mm jaw, 25cm



ヤコブソンヘビータ입持針器

Round handle, Tungsten carbide coated tips, Straight, Multilock

1802315M 2.2 × 18mm jaw, 20cm

1802317M 2.2 × 18mm jaw, 23cm

1802318M 2.2 × 18mm jaw, 25cm



ヤコブソンヘビータ입持針器

Round handle, Short knurling, Tungsten carbide coated tips, Multilock, 22cm

1802264M Straight, 0.5mm tips

1802266M Straight, 0.7mm tips

1802268M Straight, 1.0mm tips

1802270M Straight, 1.5mm tips

1802274M Curved, 0.5mm tips

1802276M Curved, 0.7mm tips

1802278M Curved, 1.0mm tips



マイクロ持針器


Round handle, Tungsten carbide coated tips, Straight, 21cm

1802164M 0.5mm tips, Multilock

1802164L 0.5mm tips, with lock

1802168M 1.1mm tips, Multilock

1802168L 1.1mm tips, with lock



カストロピエホマイクロ持針器

Round handle, without lock, Streamline jaws, tungsten carbide coated tips, 15cm

1802320 1.0mm tips, without lock

1802320L 1.0mm tips, with lock



ヤコブソンマイクロ持針器

Round handle, Streamline jaws, Tungsten carbide coated tips

1802321 1.0mm tips, 18.5cm, without lock

1802321L 1.0mm tips, 18.5cm, with lock

1802322 1.0mm tips, 20cm, without lock

1802322L 1.0mm tips, 20cm, with lock

1802323 1.1mm tips, 21cm, without lock

1802323L 1.1mm tips, 21cm, with lock

1802324 1.1mm tips, 23cm, without lock

1802324L 1.1mm tips, 23cm, with lock

1802325 1.1mm tips, 25cm, without lock

1802325L 1.1mm tips, 25cm, with lock

ヤコブソン持針器

(For large needles), Round handle, Streamline jaws, Tungsten carbide coated tips

1802330 1.5mm tips, 19cm, without lock

1802330L 1.5mm tips, 19cm, with lock

1802332 1.5mm tips, 23cm, without lock

1802332L 1.5mm tips, 23cm, with lock

1802333 1.5mm tips, 25cm, without lock

1802333L 1.5mm tips, 25cm, with lock



1802332L

1802322L

1802325L

マイクロ持針器

Flat handle, Tungsten carbide coated tips, Straight

- 1802009** Without lock, 0.5×8mm tips, 16cm
- 1802010** With lock, 0.5×8mm tips, 16cm
- 1802011** Without lock, 0.5×9mm tips, 18.4cm
- 1802012** With lock, 0.5×9mm tips, 18.4cm



マイクロ持針器

Flat handle, Tungsten carbide coated tips, Straight

- 1802013** Without lock, 0.7mm tips, 18.5cm
- 1802013L** With lock, 0.7mm tips, 18.5cm
- 1802015** Without lock, 1.0mm tips, 18.5cm
- 1802015L** With lock, 1.0mm tips, 18.5cm
- 1802017** Without lock, 1.5mm tips, 18.5cm
- 1802017L** With lock, 1.5mm tips, 18.5cm

- 1802018** Without lock, 0.5mm tips, 20cm
- 1802018L** With lock, 0.5mm tips, 20cm
- 1802019** Without lock, 0.7mm tips, 20cm
- 1802019L** With lock, 0.7mm tips, 20cm
- 1802020** Without lock, 1.0mm tips, 20cm
- 1802020L** With lock, 1.0mm tips, 20cm
- 1802021** Without lock, 1.5mm tips, 20cm
- 1802021L** With lock, 1.5mm tips, 20cm

- 1802022** Without lock, 0.5mm tips, 21cm
- 1802022L** With lock, 0.5mm tips, 21cm
- 1802023** Without lock, 0.7mm tips, 21cm
- 1802023L** With lock, 0.7mm tips, 21cm
- 1802024** Without lock, 1.0mm tips, 21cm
- 1802024L** With lock, 1.0mm tips, 21cm
- 1802025** Without lock, 1.5mm tips, 21cm
- 1802025L** With lock, 1.5mm tips, 21cm

- 1802026** Without lock, 0.5mm tips, 24cm
- 1802026L** With lock, 0.5mm tips, 24cm
- 1802027** Without lock, 0.7mm tips, 24cm
- 1802027L** With lock, 0.7mm tips, 24cm
- 1802028** Without lock, 1.0mm tips, 24cm
- 1802028L** With lock, 1.0mm tips, 24cm
- 1802029** Without lock, 1.5mm tips, 24cm
- 1802029L** With lock, 1.5mm tips, 24cm



ヤサーギルマイクロ持針器

Bayonet style, Flat handle, Tungsten carbide coated tips,
0.5mm tips,
Working distance: 78mm
Length: 18.5cm

- 1802401** Straight, without lock
- 1802402** Straight, with lock
- 1802403** Curved, without lock
- 1802404** Curved, with lock

Bayonet style, Flat handle, Tungsten carbide coated tips,
0.5mm tips,
Working distance: 85mm
Length: 20cm

- 1802405** Straight, without lock
- 1802406** Straight, with lock
- 1802407** Curved, without lock
- 1802408** Curved, with lock

Bayonet style, Flat handle, Tungsten carbide coated tips,
0.5mm tips,
Working distance: 95mm
Length: 22.5cm


- 1802409** Straight, without lock
- 1802410** Straight, with lock
- 1802411** Curved, without lock
- 1802412** Curved, with lock



ロートン持針器

Bayonet style, Tungsten carbide coated tips, Round handle

- 1802501** 0.9mm tips, Working distance 89mm, without lock, Straight, 21.6cm
- 1802503** 0.9mm tips, Working distance 89mm, without lock, Upwards curved, 21.6cm
- 1802505** 0.9mm tips, Working distance 95mm, without lock, Straight, 23cm
- 1802507** 0.9mm tips, Working distance 95mm, without lock, Upwards curved tip, 23cm
- 1802506** 0.9mm tips, Working distance 95mm, with lock, Straight, 23cm
- 1802508** 0.9mm tips, Working distance 95mm, with lock, Upwards curved tip, 23cm



ヘガール型持針器

Ring handle, with tungsten carbide coated tips

1815012 Standard jaws, 12.5cm

1815012A Fine jaws, 12.5cm

1815012B Strong jaws, 12.5cm

1815013 Standard jaws, 13cm

1815013A Fine jaws, 13cm

1815013B Strong jaws, 13cm

1815014 Standard jaws, 14cm

1815014A Fine jaws, 14cm

1815014B Strong jaws, 14cm

1815016 Standard jaws, 16cm


1815016A Fine jaws, 16cm

1815016B Strong jaws, 16cm

1815018 Standard jaws, 18cm

1815018A Fine jaws, 18cm

1815018B Strong jaws, 18cm




ライダー持針器

1800730 1.4 ×10mm tungsten carbide coated tips, 15cm

1800732 1.4 ×11mm tungsten carbide coated tips, 17cm

1800734 1.4 ×12mm tungsten carbide coated tips, 19cm




マイクロ持針器

1800735 Extra fine tungsten carbide coated tips, 19cm

フェル持針器

1800737 2mm tungsten carbide inlay jaws, 21cm

1800740 2mm tungsten carbide inlay jaws, 30cm



1800735



1800740

クライル・ウッド持針器

Tungsten carbide inlay jaws

1800762 Jaw 20 × 2mm, 20cm

1800763 Jaw 20 × 3mm, 20cm

1800764 Jaw 20 × 4mm, 20cm

1800760 Jaw 21 × 2mm, 23cm

胸骨針用持針器


1800766 Tungsten carbide inlay jaws,
Jaw 20 × 6mm, 20cm



ドバーキ ダブルアクション持針器
Tungsten carbide inlay jaws

1800824 Jaw 20.7 × 2.2mm, 24cm

1800828 Jaw 20.7 × 2.5mm, 28cm



1800828

オルセン・ヘガール剪刃付持針器

With tungsten carbide coated tips, Combined with suture scissors

1816012 Standard jaws, 12.5cm

1816012A Fine jaws, 12.5cm

1816013 Standard jaws, 13cm

1816013A Fine jaws, 13cm

1816014 Standard jaws, 14cm

1816014A Fine jaws, 14cm

1816016 Standard jaws, 16cm

1816016A Fine jaws, 16cm

1816018 Standard jaws, 18cm

1816018A Fine jaws, 18cm



1816012A



1816012

外科剪刀

- 1812113** Straight, Sharp-sharp, 14cm
- 1812133** Straight, Sharp-blunt, 14cm
- 1812153** Straight, Blunt-blunt, 14cm
- 1812213** Curved, Sharp-sharp, 14cm
- 1812233** Curved, Sharp-blunt, 14cm
- 1812253** Curved, Blunt-blunt, 14cm

Sharp - Sharp




Sharp - Blunt



Blunt - Blunt



1812153



リスター包帯剪刀

- 1812110** 4cm long blades, 15cm long



ミニヴァンナス剪刀

1803001 Sharp pointed tips, 7mm from pivot to tip
Cutting edge 3mm, 7.5cm



ヴァンナス剪刀

Sharp pointed tips, Cutting edge 3mm

1803004 Straight, 7.8cm

1803005 Curved, 7.8cm

1803006 Angled, 7.8cm




1803004

ノイス剪刀


Sharp points, 21mm from pivot to tip, 11.4cm

1803007 Curved

1803008 Straight



1803007



1803008

ウエストコット剪刀

Blunt tips, 21mm from pivot to tip, 115mm long

1803009 Curved

1803010 Straight




1803009

ヤコブソンマイクロ剪刀

Flat handle, 12mm blades, 17.8cm



- 1803015** Straight
- 1803016** Curved
- 1803017** Angled on flat 25° blades
- 1803018** Angled on flat 45° blades
- 1803019** Angled on flat 60° blades
- 1803020** Angled on flat 90° blades



マイクロ剪刀

Round handle, 13mm blades, 14cm

- 1803041** Curved delicate (0.3mm tips)
- 1803042** Straight delicate (0.3mm tips)
- 1803043** Curved (0.7mm tips)
- 1803044** Straight (0.7mm tips)



マイクロ剪刀

Round handle, 15mm blades

- 1803081** Curved delicate (0.3mm tips), 18cm
- 1803082** Straight delicate (0.3mm tips), 18cm
- 1803083** Curved (0.7mm tips), 18cm
- 1803084** Straight (0.7mm tips), 18cm

- 1803111** Curved delicate (0.4mm tips), 21cm
- 1803112** Straight delicate (0.4mm tips), 21cm
- 1803113** Curved (0.7mm tips), 21cm
- 1803114** Straight (0.7mm tips), 21cm

- 1803101** Curved delicate (0.4mm tips), 20cm
- 1803102** Straight delicate (0.4mm tips), 20cm
- 1803103** Curved (0.7mm tips), 20cm
- 1803104** Straight (0.7mm tips), 20cm


- 1803131** Curved delicate (0.4mm tips), 23cm
- 1803132** Straight delicate (0.4mm tips), 23cm
- 1803133** Curved (0.7mm tips), 23cm
- 1803134** Straight (0.7mm tips), 23cm



ヤコブソン・マイクロ剪刀

Flat handle, Fine blades, 16.5cm

- 1803241** 25° angle
- 1803242** 45° angle
- 1803243** 60° angle
- 1803244** 90° angle
- 1803245** 125° angle
- 1803246** V-neck, 125° angle



ポッツ・マイクロ剪刀

Flat handle, Short fine blades

1803021 25° angle, 18cm

1803022 45° angle, 18cm

1803023 60° angle, 18cm

1803024 90° angle, 17.8cm

1803025 125° angle, 17.8cm

1803026 25° angle, 20.3cm

1803027 45° angle, 20.3cm

1803028 60° angle, 20.3cm

1803029 90° angle, 20cm

1803030 125° angle, 20cm



1803026



1803027



1803028



1803029



1803030

ポッツ・マイクロ剪刀

Round handle, Short fine blades

1803201	25° angle, 17cm	1803211	25° angle, 20cm	1803216	25° angled, 18cm
1803202	45° angle, 17cm	1803212	45° angle, 20cm	1803217	45° angled, 18cm
1803203	60° angle, 17cm	1803213	60° angle, 20cm	1803218	60° angled, 18cm
1803204	90° angle, 16.5cm	1803214	90° angle, 19.5cm	1803219	90° angled, 17.8cm
1803205	125° angle, 16.5cm	1803215	125° angle, 19.5cm	1803220	125° angled, 17.8cm



ヤコブソン・マイクロ剪刀

Round handle, 1 blade with bead tip, Short fine blades

- 1803221** 25° angle, 18cm
- 1803222** 45° angle, 18cm
- 1803223** 60° angle, 18cm
- 1803224** 90° angle, 17.8cm
- 1803225** 125° angle, 17.8cm



ポッツ・マイクロ剪刀


Bayonet flat handle, 12mm blades

- 1803302** 45° angle, 19cm
- 1803305** 125° angle, 18.5cm
- 1803310** 125° reverse angle cutting blades (left handed), 18.5cm
- 1803312** 45° angle, 1 blade with bead tip, 19cm
- 1803315** 125° angle, 1 blade with bead tip, 18.5cm



ヤコブソン・ヴァンナスマイクロ剪刀

- 1803031** 17cm



ヤサーギル・マイクロ剪刀

Bayonet flat handle

1803341 Straight, 13mm blades, 18.5cm

1803341S Straight, Serrated upper blade,
13mm blades, 18.5cm

1803342 Curved, 13mm blades, 18.5cm

1803342S Curved, Serrated upper blade,
13mm blades, 18.5cm

1803343 Straight, 14mm blades, 20cm

1803343S Straight, Serrated upper blade,
14mm blades, 20cm

1803344 Curved, 14mm blades, 20cm

1803344S Curved, Serrated upper blade,
14mm blades, 20cm

1803345 Straight, 15mm blades, 23cm


1803345S Straight, Serrated upper blade,
15mm blades, 23cm

1803346 Upwards Curved, 15mm blades, 23cm

1803346S Upwards Curved, Serrated upper blade,
15mm blades, 23cm

1803347 Downwards Curved, 15mm blades, 23cm

1803347S Downwards Curved, Serrated upper blade,
15mm blades, 23cm



マリス・バイオネットマイクロ剪刀

Flat handle, 23cm

1803351 Straight, Regular tips

1803351S Straight, Regular tips,
Serrated upper blade

1803352 Straight, Fine tips

1803352S Straight, Fine tips,
Serrated upper blade

1803353 Straight, Fine tenotomy tips

1803353S Straight, Fine tenotomy tips,
Serrated upper blade

1803354 Curved up, Regular tips

1803354S Curved up, Regular tips,
Serrated upper blade


1803355 Curved up, Fine tips

1803355S Curved up, Fine tips,
Serrated upper blade

1803356 Curved down, Regular tips

1803356S Curved down, Regular tips,
Serrated upper blade

1803351



1803352S

バイオネットマイクロ剪刀

1803361 Straight, Regular tips, 25.7cm

1803361S Straight, Regular tips,
Serrated upper blade, 25.7cm

1803362 Curved, Regular tips, 25.7cm

1803362S Curved, Regular tips,
Serrated upper blade, 25.7cm

1803361



バイオネットマイクロ剪刀

Flat handle, Bayonet style


- 1803371** 5mm blades, Straight, 18.5cm
- 1803372** 5mm blades, Upwards, curved, 18.5cm
- 1803381** 5mm blades, Straight, 22.5cm
- 1803382** 5mm blades, Upwards curved, 22.5cm



マリス・バイオネットマイクロ剪刀

Round handle, 23cm


- 1803500** Straight, Regular tips
- 1803500S** Straight, Regular tips,
Serrated upper blade
- 1803501** Curved up, Regular tips
- 1803501S** Curved up, Regular tips,
Serrated upper blade
- 1803502** Straight, Fine tenotomy tips
- 1803502S** Straight, Fine tenotomy tips,
Serrated upper blade



レザノ剪刀 三連関節式
With 3 articulations

1800924 24cm

1800928 28cm





1800924



1800928


外科剪刀

- 1800930 Plastic-Jabaley Scissors, Straight, 12cm
- 1800932 Plastic-Jabaley Scissors, Curved, 12cm
- 1800942 Baby Metzenbaum Scissors, Straight, 12.5cm
- 1800944 Baby Metzenbaum Scissors, Curved, 12.5cm
- 1800947 Metzenbaum Scissors, Straight, 18cm
- 1800948 Metzenbaum Scissors, Curved, 18cm
- 1800949 Metzenbaum Scissors, Straight, 20cm
- 1800950 Metzenbaum Scissors, Curved, 20cm
- 1800956 Metzenbaum Scissors, Curved, 30.5cm



メーヨー剪刀

- 1800962 Mayo Harrington Scissors**
Straight, 20cm
- 1800963 Mayo Harrington Scissors**
Straight, 23cm
- 1800966 Mayo Harrington Scissors**
Curved, 20cm
- 1800972 Mayo Harrington Scissors**
Curved, 30.5cm



1800962

ヤサーギル・マイクロ結紮ガイド

- 1808401 Slightly curved, 18.5cm
- 1808402 Curved, 18.5cm
- 1808403 90° angled, 18.5cm
- 1808404 135° angled, 18.5cm



- 1808408 ヤサーギル・マイクロ血管メス 18.5cm



ヤサーギル・マイクロ剥離子

- 1808409 Angled, 18.5cm
- 1808410 Round tip, 18.5cm
- 1808411 Slightly curved, 18.5cm
- 1808412 Curved, 18.5cm



- 1808413 ヤサーギルマイクロディセクター 18.5cm



- 1808414 ヤコブソン・マイクロプローブ 18.5cm



1808416 ヤコブソン・マイクロ結紮フォーク ,18.5cm



クライエンプル血管・神経鉤

1808418 1.0mm diameter, Hook depth 3mm, Probe pointed small,18.5cm

1808419 1.0mm diameter, Hook depth 3.8mm, Probe pointed large,18.5cm

1808420 Without probe pointed small, Hook depth 2.5mm,18.5cm

1808421 Without probe pointed large, Hook depth 4mm,18.5cm



1808422 マイクロフック 18cm



1808431 3mm diameter, Straight, 21cm



キュレット



バイオネットマイクロ血管メス

1808700 Bayonet style, Sharp tip,18.5cm

1808701 Bayonet style, Blunt tip,18.5cm



ヤサーギルマイクロくも膜メス

- 1808702** Upwards curved, 18.5cm
- 1808703** Downwards curved, 18.5cm
- 1808712** Upwards curved, 21.5cm
- 1808713** Downwards curved, 21.5cm



1808702



1808703

1808715 マイクロ鋭匙

22cm



ヤサーギルマイクロ剥離子

Bayonet style

- 1808704** Downwards curved, 18.5cm
- 1808705** Upwards curved, 18.5cm
- 1808706** Slight curved, 18.5cm



1808704





1808705




1808706

バイオネットディセクター


- 1808716** 1.15mm wide tip, Small, 23cm
1808717 1.7mm wide tip, Medium, 23cm
1808718 2.1mm wide tip, Large, 23cm



1808716




1808717



1808718

バイオネットエレベーター

- 1808719** Oval tip, 2.75mm wide, 22.5cm




バイオネット神経鉤

- 1808720** Bulbous tip, 23cm
1808721 Angle 90°, 0.5mm tip, 23cm



1808720



1808720



1808721

バイオネットマイクロリングキュレット

- 1808722** Angled up 90°, 3.0mm diameter, 25.5cm
- 1808723** Angled down 90°, 3.0mm diameter, 25.5cm
- 1808724** Angled up 45°, 2.25mm diameter, 25.5cm
- 1808725** Angled up 90°, 3.0mm diameter, 28cm
- 1808726** Angled down 40°, 3.0mm diameter, 28cm
- 1808727** Angled up 45°, 2.25mm diameter, 28cm



ギャレット血管拡張器

Malleable shaft

- 1808748A** 0.5mm, 14cm
- 1808748B** 1.25mm, 14cm
- 1808748C** 1.75mm, 14cm
- 1808748D** 2.0mm, 14cm

- 1808749A** 0.5mm, 21cm
- 1808749B** 1.25mm, 21cm
- 1808749C** 1.75mm, 21cm
- 1808749D** 2.0mm, 21cm



カルボンティ血管鉗

1808845 Angle 45°

1808890 Angle 90°

1808925 Angle 125°

1808845



1808890



1808925



I MA用心外膜レトラクター

1874047 47mm

1874063 54mm




1822200 ウッド両頭鉤

Blade A: 28mm D × 20mm W,


Blade B: 40mm D × 20mm W,

Overall length 315mm




クーリー心房鉤（僧帽弁用鉤）

- 1822236** Large, Platypus, 36mm blade width, 27cm
1822221 Medium, Platypus, 23mm blade width, 25.5cm
1822231 Small, Platypus, 23mm blade width, 23.5cm



- 1874043** 新生児用胸骨用開創器
With interchangeable blades (2 pair), Arm length 52mm




- 1874083** フィノシェット胸骨用開創器 新生児
Arm length 52mm



1874065 胸骨開胸器
Blades 12mm D × 25mm W,
Arm length 83mm



1874028 フィノシェット肋骨開胸器
Blades 30mm D × 32mm W,
Arm length 83mm



モース胸骨開胸器

1874066 Child size


Total opening	160mm
Blade width	19mm
Blade depth	19mm
Arm length	156mm

1874067 Adult size


Total opening	200mm
Blade width	28mm
Blade depth	25mm
Arm length	200mm

1874090 Adult size


Total opening	200mm
Blade width	28mm
Blade depth	38mm
Arm length	200mm



1874066



1874067



1874090

1874170 デュボスト開胸器


Arm length 185mm, Curved, with interchangeable blades (3 pair)

Consists of:



74170-30 Small blade

74170-40 Medium blade



74170-50 Medium blade





74170-30



74170-40



74170-50




1874211 カルボンティ開胸器
Arm length is 200mm, Curved

Consists of:

74211-15 15mm malleable narrow blade

74211-25 25mm malleable wide blade


74255-00 Median blade with slider



74211-15



74211-25



74255-00

1874252 IMA開胸器 内胸動脈用
Arm length 170mm

Consists of:

- 74252-20** 20mm blades
- 74252-25** 25mm blades
- 74252-28** 28mm hooks
- 74252-32** 32mm hooks
- 74252-38** 38mm hooks

